



## PRODUCT SPECIFICATION

catalog number:	<b>A46291</b>
product description:	HE4 Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Ovarian carcinoma cells
source:	Murine
form:	Liquid
purity:	> 90%
buffer:	10mM PBS, pH 7.4
preservative:	0.1% P300
concentration:	[lot specific] mg/mL
manufacture date:	[lot specific]
expiration date:	[lot specific]
storage:	2-8°C short term (two weeks); -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	There is no known hazard/biohazard associated with this product.