



PRODUCT SPECIFICATION

catalog number:	A46400
product description:	CA 19.9 Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Purified CA 19.9
source:	Culture supernatant
form:	Liquid
purification method:	Ion exchange chromatography
buffer:	15 mM phosphate, 0.15 M NaCl, pH 7.2
preservative:	0.09% NaN ₃
isotype:	IgA
concentration:	[lot specific] mg/mL (based on A ^{1%} _{280nm} = 14.0)
specificity:	This monoclonal antibody is specific for human CA 19.9.
pairing:	This antibody can be used as both the capture and detection antibody in sandwich EIA.
storage:	2-8°C short term; -70°C long term (up to 12 months) avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.