



## PRODUCT SPECIFICATION

catalog number:	<b>A46590</b>								
product description:	Dengue Virus NS1 Monoclonal Antibody, Purified								
lot number:	[lot specific]								
immunogen:	Recombinant dengue virus, type NS1								
source:	Murine ascites								
form:	Liquid								
purification method:	DEAE chromatography								
purity:	>95% by SDS-PAGE and HPLC								
buffer:	0.015M potassium phosphate buffer, 0.85% NaCl, pH 7.2								
preservative:	0.05% NaN <sub>3</sub>								
isotype:	IgG <sub>1</sub> , kappa								
concentration:	[lot specific] mg/mL (based on A <sup>1%</sup> <sub>280nm</sub> = 14.0)								
specificity:	This antibody is specific for dengue virus NS1, type 2 & 3.								
cross-reactivity:	<table><tr><td>▪ Type 1, NS1</td><td>0.0%</td><td>▪ Type 3, NS1</td><td>150%</td></tr><tr><td>▪ Type 2, NS1</td><td>100%</td><td>▪ Type 4, NS1</td><td>0.0%</td></tr></table>	▪ Type 1, NS1	0.0%	▪ Type 3, NS1	150%	▪ Type 2, NS1	100%	▪ Type 4, NS1	0.0%
▪ Type 1, NS1	0.0%	▪ Type 3, NS1	150%						
▪ Type 2, NS1	100%	▪ Type 4, NS1	0.0%						
pairing:	This antibody can be used in sandwich ELISA in conjunction with A46582.								
storage:	2-8°C short term (2 years); -20°C long term avoiding multiple freeze thaw cycles.								
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.								