



PRODUCT SPECIFICATION

catalog number:	A47013
product description:	RSV Monoclonal Antibody
lot number:	[lot specific]
immunogen:	RSV strain; long
source:	Murine ascites
form:	Liquid
purification method:	Molecular exclusion chromatography
purity:	[lot specific]
buffer:	Borate saline buffer, pH 9.2
preservative:	0.09% NaN ₃
isotype:	IgM
concentration:	[lot specific] mg/mL (A _{280nm})
specificity:	RSV
expiration date:	[lot specific]
storage:	2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.