



## PRODUCT SPECIFICATION

**CATALOG NUMBER: A47030**

PRODUCT:	Respiratory syncytial virus (RSV) Monoclonal Antibody
IMMUNOGEN:	Purified RSV virus, Strain Long
PURIFICATION:	Chromatography on protein G Sepharose.
BUFFER:	PBS, pH 7.4
PRESERVATIVE:	0.1% sodium azide
SUBCLASS:	IgG <sub>2b</sub>
STORAGE CONDITIONS:	2 - 8°C
SPECIFICITY:	This antibody is specific To F protein
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.