



PRODUCT SPECIFICATION

catalog number:	A48061
product description:	Strep A Monoclonal Antibody, Purified
lot number:	[lot specific]
source:	Rabbit; culture medium
form:	Liquid
purification method:	Protein A chromatography
purity:	>90%
buffer:	0.01 M phosphate buffered saline, pH 7.2
preservative:	0.1% NaN ₃
isotype:	IgG ₁
concentration:	[lot specific] mg/mL
specificity:	Specific to the Group A Strep carbohydrate.
cross-reactivity:	No cross reactivity with: Streptococci groups B, C, D, F & G + other bacteria
pairing:	This antibody is self-pairing and can also be paired with G47091041 or G47090041 in lateral flow.
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.