



PRODUCT SPECIFICATION

catalog number:	A480720
product description:	Staphylococcal Enterotoxin B Monoclonal Antibody
lot number:	[lot specific]
source:	Murine
form:	Liquid
purification method:	Protein A chromatography
purity:	>90 %
buffer:	0.01 M Phosphate buffered saline, pH 7.2
preservative:	0.01% NaN ₃
concentration:	[lot specific] mg/ml
specificity:	Specific to enterotoxin B of Staphylococcal aureus
pairing:	This antibody can be used in sandwich ELISA as the capture antibody in conjunction with A48073 as the detector.
storage:	2-8°C; DO NOT FREEZE
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.