



## PRODUCT SPECIFICATION

catalog number:	<b>A48120</b>
product description:	Chikungunya Virus (CHKV) Monoclonal Antibody, Purified
lot number:	[lot specific]
source:	Murine
form:	Liquid
purification method:	Protein A chromatography
buffer:	0.01M PBS, pH7.2
preservative:	0.1% NaN <sub>3</sub>
isotype:	IgG <sub>1</sub>
concentration:	[lot specific] mg/mL
cross-reactivity:	This antibody shows no cross-reactivity with other alphaviruses.
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.