



## PRODUCT SPECIFICATION

catalog number:	<b>A48127</b>
product description:	Chikungunya Virus (CHKV) Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	BEI inactivated vaccine strain of CHKV (181/25)
source:	Tissue culture supernatant
form:	Liquid
purification method:	Protein A chromatography
buffer:	0.15M PBS
preservative:	none
isotype:	IgG <sub>1</sub> , kappa
concentration:	[lot specific] mg/mL
cross-reactivity:	<ul style="list-style-type: none"><li>▪ RR64 virus 0%</li><li>▪ O'nyong nyong virus 0%</li><li>▪ Mayaro virus 0%</li></ul>
storage:	-20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	There is no known hazard or biohazard associated with this product.