



## PRODUCT SPECIFICATION

**CATALOG NUMBER: A49010**

PRODUCT:	Thyroxine binding globulin (TBG) Monoclonal Antibody
FORM:	Protein A purified immunoglobulin from mouse ascites
PURITY:	>90% by SDS-PAGE
BUFFER:	10 mM sodium phosphate buffer, 0.154M NaCl, pH 7.4
PRESERVATIVE:	None
AFFINITY CONSTANT:	$1 \times 10^9$
SUBCLASS:	IgG <sub>2a</sub> kappa
STORAGE CONDITIONS:	< -10°C
REFERENCE ASSAY:	RIA
COUNTRY OF ORIGIN:	U.S.A.
HAZARD/BIOHAZARD:	There is no known biohazard associated with this product.