



## PRODUCT SPECIFICATION

**CATALOG NUMBER: A49400**

PRODUCT:	Cystatin C Monoclonal Antibody
LOT NUMBER:	Reference
PURIFICATION:	Purified immunoglobulin from BALB/C mouse ascites (Affinity chromatography)
BUFFER:	15mM Phosphate buffer containing 0.15M NaCl, pH 7.2
PRESERVATIVE:	0.1% sodium azide
SUBCLASS:	IgG <sub>1k</sub>
CONCENTRATION:	Reference (by $A_{280nm}^{1\%} = 14.0$ )
STORAGE CONDITIONS:	2-8°, Maintain at -70°C long term. Avoid repeated freeze-thaw cycles.
APPLICATIONS:	Sandwich ELISA; can be used as capture antibody when paired with A49410.
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precaution. There is no known biohazard associated with this product.