



## PRODUCT SPECIFICATION

catalog number:	<b>A49410</b>
product description:	Cystatin C Monoclonal Antibody, Purified
lot number:	[lot specific]
source:	Murine ascites
form:	Colorless liquid
purification method:	Affinity chromatography
buffer:	15mM phosphate buffer, 0.15M NaCl, pH 7.2
preservative:	0.09% NaN <sub>3</sub>
isotype:	IgG <sub>2a</sub> , kappa
concentration:	[lot specific] mg/mL (based on A <sup>1%</sup> <sub>280nm</sub> = 14.0)
specificity:	This antibody is specific for human cystatin C.
application/notes:	This monoclonal can be used in a sandwich ELISA as the conjugate antibody in combination with A49400 (as the capture antibody).
expiration date:	[lot specific]
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.