



PRODUCT SPECIFICATION

catalog number:	A50091
product description:	K2 Monoclonal Antibody, Purified
lot number:	[lot specific]
form:	Liquid
buffer:	PBS, pH 7.4
preservative:	0.02% NaN ₃
concentration:	[lot specific] mg/mL
pairing:	This antibody pairs well with the K2 antigen, J50091.
expiration date:	[lot specific]
storage:	2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.