



## PRODUCT SPECIFICATION

**CATALOG NUMBER: A51020**

PRODUCT:	Tobramycin Monoclonal Antibody
IMMUNOGEN:	Tobramycin-BSA
FORM:	Purified IgG
PURIFICATION:	Protein A
BUFFER:	10 Mm Phosphate, 150 mM NaCl, pH 7.4
PRESERVATIVE:	0.1% Sodium Azide
AFFINITY CONSTANT:	$1 \times 10^{10}$
SUBCLASS:	IgG <sub>1</sub> ,Kappa
STORAGE CONDITIONS:	2-8°C
HAZARD/BIOHAZARD:	This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.