



PRODUCT SPECIFICATION

CATALOG NUMBER: A51950

PRODUCT:	Benzodiazepine Monoclonal Antibody
IMMUNOGEN:	BZO-BSA Conjugate
PURITY:	>95% by SDS Page
PURIFICATION:	Purified IgG from Mouse Ascites
BUFFER:	Phosphate Buffer Saline, pH 7.0
PRESERVATIVE:	0.02% sodium azide
SUBCLASS:	IgG _{1K}
STORAGE CONDITIONS:	2-8°C short term -20°C long term
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precaution. There is no known biohazard associated with this product.