



PRODUCT SPECIFICATION

catalog number:	A52015
product description:	Benzoylcegonine Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Benzoylcegonine-BTG
source:	Murine
form:	Liquid
purification method:	Protein A chromatography
purity:	≥95% by SDS-PAGE
buffer:	PBS, pH 7.4
preservative:	0.1% NaN ₃
concentration:	[lot specific] mg/mL
expiration date:	Three years from date of receipt
storage:	2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.
country of origin:	USA
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.