



PRODUCT SPECIFICATION

catalog number:	A52076
product description:	THC Monoclonal Antibody
Immunogen	THC-BSA
lot number:	[lot specific]
source:	Murine
Purification	DEAE chromatography
buffer:	0.015 M Potassium Phosphate Buffer, 0.85% NaCl, pH 7.2
preservative:	0.05% NaN ₃
isotype:	IgG ₁ Kappa
concentration:	[lot specific] mg/mL
specificity:	THC and Metabolites
pairing:	This antibody can be used in sandwich test as the detection antibody in conjunction with V50000527 (as the capture).
notes/applications:	Lateral Flow, EIA
storage:	2-8°C short term (2 years); -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.