



PRODUCT SPECIFICATION

catalog number:	A52140
product description:	Phencyclidine (PCP) Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	PCP-BSA
source:	Murine
form:	Liquid
buffer:	0.05M Phosphate, 0.85% NaCl, pH 7.15
preservative:	0.05% NaN ₃
isotype:	IgG ₁
concentration:	[lot specific] mg/mL
specificity:	This monoclonal antibody is specific for phencyclidine.
cross-reactivity:	There is no apparent cross reactivity to benzoylcegonine, opiates, cannabinoids or amphetamines.
expiration date:	[lot specific]
storage:	2-8°C (≤ two weeks); -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.