



PRODUCT SPECIFICATION

catalog number:	A52148503P
product description:	Phencyclidine (PCP) Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	PCP-BTG
source:	Murine
form:	Liquid
purification method:	DEAE chromatography
purity:	>95% by HPLC and SDS-PAGE
buffer:	0.015M potassium phosphate buffer, 0.85% NaCl, pH 7.2
preservative:	0.05% NaN ₃
isotype:	IgG ₁
concentration:	[lot specific] mg/mL (based on A ^{1%} _{280nm} = 14.0)
specificity:	This monoclonal antibody is specific for phencyclidine.
applications:	Useful for quantitative PCP assays by EIA
storage:	2-8°C short term (1year); -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.