



PRODUCT SPECIFICATION

catalog number:	A52210
product description:	Buprenorphine (BUP) Monoclonal Antibody, Purified
lot number:	[lot specific]
source:	Murine
form:	Liquid
buffer:	15mM phosphate, 0.85% NaCl, pH7.2
preservative:	0.1% NaN ₃
isotype:	IgG
concentration:	[lot specific] mg/mL
storage:	2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.