



**PRODUCT SPECIFICATION**

**CATALOG NUMBER: A52310**

PRODUCT:	Ketamine Monoclonal Antibody
IMMUNOGEN:	Ketamine-BgG
BUFFER:	0.05M Phosphate 0.85% Sodium Chloride, pH 7.15
PRESERVATIVE:	0.05% NaN <sub>3</sub>
SUBCLASS:	IgG <sub>1</sub>
SPECIFICITY:	Ketamine
STORAGE CONDITIONS:	2-8°C Short Term, Long Term - Frozen