



## PRODUCT SPECIFICATION

catalog number:	<b>A52520</b>
product description:	Lysergic acid diethylamide (LSD) Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	LSD-BgG
source:	Murine
form:	Liquid
buffer:	0.05M K <sub>3</sub> PO <sub>4</sub> , 0.85% NaCl, pH 7.5
preservative:	0.05% NaN <sub>3</sub>
isotype:	IgG
concentration:	[lot specific] mg/mL
storage:	2-8°C short term (≤ 2 weeks); -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.