



PRODUCT SPECIFICATION

catalog number:	A53044
product description:	Fentanyl Monoclonal Antibody
lot number:	[lot specific]
immunogen:	Fentanyl-BSA
source:	Murine ascites
form:	Liquid
purification method:	Protein A affinity chromatography
purity:	>95%
buffer:	50 mM NaCl, 10 mM PBS, pH 8.0
preservative:	0.1% NaN ₃
isotype:	IgG ₁
concentration:	[lot specific] mg/mL
pairing:	This antibody when used as the label forms a match pair with V53044 (as the coating).
notes/applications:	ELISA or Lateral flow
storage:	2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.