



PRODUCT SPECIFICATION

| | |
|----------------------|---|
| catalog number: | A53136 |
| product description: | Amphetamine Monoclonal Antibody, Purified |
| lot number: | [lot specific] |
| immunogen: | Amp-BSA conjugate |
| source: | Murine ascites |
| form: | Liquid |
| purity: | >95% by SDS-PAGE |
| buffer: | PBS, pH 7.0 |
| preservative: | 0.02% NaN ₃ |
| isotype: | IgG ₁ |
| concentration: | [lot specific] mg/mL (based on A280nm, 1mg/ml, 1cm = 1.4) |
| storage: | 2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles. |
| hazard/biohazard: | This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product. |