



PRODUCT SPECIFICATIONS

CATALOG NUMBER: A53149

PRODUCT:	EDDP Monoclonal Antibody
IMMUNOGEN:	EDDP-BSA Conjugate
FORM:	Purified IgG from mouse ascites
PURITY:	>95% by SDS-Page
BUFFER:	Phosphate buffer saline, pH 7.0
PRESERVATIVE:	0.02% sodium azide
SUBCLASS:	IgG ₁
STORAGE CONDITIONS:	Short term: 2-8°C; Long term: -20 °C
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.