



## PRODUCT SPECIFICATION

catalog number:	<b>A53330</b>
product description:	MDMA Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	MDMA-BTG
source:	Murine
form:	Liquid
purification method:	DEAE chromatography
purity:	>95% by HPLC & SDS-PAGE
buffer:	0.015M potassium phosphate, 0.85% NaCl, pH 7.2
preservative:	0.05% NaN <sub>3</sub>
isotype:	IgG <sub>1, kappa</sub>
concentration:	[lot specific] mg/mL (based on A <sup>1%</sup> <sub>280nm</sub> = 14.0)
specificity:	This antibody is specific for MDMA.
cross-reactivity:	▪ Amphetamine <0.01%      ▪ Methamphetamine 0.01%
storage:	2-8°C (two years): -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.