



## PRODUCT SPECIFICATION

catalog number:	<b>A53370574P</b>
product description:	MDMA Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	MDMA-BSA conjugate
source:	Murine ascites
form:	Liquid
purity:	95% by affinity chromatography
buffer:	50mM NaCl and 10mM Gly-Hcl
preservative:	0.05% NaN <sub>3</sub>
isotype:	IgG <sub>1</sub>
concentration:	[lot specific] mg/mL
specificity:	This antibody is specific for MDMA.
expiration date:	[lot specific]
storage:	2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.