



PRODUCT SPECIFICATIONS

CATALOG NUMBER: A53470

PRODUCT:	Oxycodone Monoclonal Antibody
IMMUNOGEN:	Oxycodone-BSA Conjugate
FORM:	Purified IgG from mouse ascites
PURITY:	>95% by SDS-Page
BUFFER:	Phosphate Buffer Saline, pH 7.0
PRESERVATIVE:	0.02% Sodium Azide
SUBCLASS:	IgG ₁
STORAGE CONDITIONS:	Short term: 2-8°C Long term: -20°C
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precaution. There is no known biohazard associated with this product.