



## PRODUCT SPECIFICATION

catalog number:	<b>A53480</b>
product description:	Oxycodone Monoclonal Antibody, Purified
lot number:	[lot specific]
source:	Murine
form:	Liquid
purification method:	Protein A chromatography
buffer:	PBS, pH 7.4
preservative:	0.1 % NaN <sub>3</sub>
concentration:	[lot specific] mg/mL
storage:	2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.