



## PRODUCT SPECIFICATION

**CATALOG NUMBER: A53490**

PRODUCT:	Oxycodone Monoclonal Antibody
IMMUNOGEN:	Oxycodone-BSA Conjugate
FORM:	Purified immunoglobulin from mouse ascites
PURITY:	>90%. Caprylic acid method and ammonium sulfate method.
BUFFER:	0.1 M PBS, pH 7.4
STORAGE CONDITIONS:	2-8°C short term; -20°C long term.
STABILITY:	Two weeks when stored at 2-8°C Two years when stored at -20°C
HAZARD/BIOHAZARD:	There is no known hazard or biohazard associated with this product.