



PRODUCT SPECIFICATION

catalog number:	A56044
product description:	Ethyl Glucuronide (EtG) Monoclonal Antibody, Purified
lot number:	[lot specific]
source:	Sheep
form:	Liquid
purity:	>90%
buffer:	PBS, pH 7.3
preservative:	0.1% NaN ₃
concentration:	[lot specific] mg/mL
pairing:	This antibody can be paired with EtG-BSA V56045 or V56046 or EtG-HRP V56048.
notes/applications:	This monoclonal is known to work well in ELISA and EMIT.
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.