



PRODUCT SPECIFICATION

CATALOG NUMBER: A57030

PRODUCT:	<u>Toxoplasma gondii</u> Monoclonal Antibody
FORM:	Liquid
PURIFICATION:	Immunoglobulin purified by chromatography on Protein G Sepharose.
BUFFER:	PBS, pH 7.4
PRESERVATIVE:	0.1% NaN ₃
SUBCLASS:	IgG2a
STORAGE CONDITIONS:	2-8°C
SPECIFICITY:	Clone reacts with p30 membrane protein in ELISA and Immunoblotting
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.