



PRODUCT SPECIFICATION

catalog number:	A58014
product description:	C-Reactive Protein (CRP) Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Purified human serum CRP
source:	Murine ascites
form:	Liquid
purity:	>98% by SDS-PAGE
buffer:	PBS, pH 7.0
preservative:	0.02% NaN ₃
isotype:	IgG _{2a}
concentration:	[lot specific] mg/mL (based on A ^{1%} _{280nm} = 14.0)
specificity:	This antibody is specific for c-reactive protein.
storage:	2 - 8°C short term (2 years); -20°C long term
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.