



## PRODUCT SPECIFICATION

catalog number:	<b>A58015</b>
product description:	CRP Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Human serum CRP
source:	Murine
form:	Liquid
purity:	>98% by SDS-PAGE
buffer:	PBS, pH 7.0
preservative:	0.02% NaN <sub>3</sub>
isotype:	IgG <sub>2a</sub>
concentration:	[lot specific] mg/mL (based on A280nm, 1 mg/ml, 1 cm = 1.40)
storage:	2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.