



PRODUCT SPECIFICATION

catalog number:	A58017
product description:	CRP Monoclonal Antibody
lot number:	[lot specific]
source:	Murine ascites
form:	Liquid
purification method:	Affinity chromatography
purity:	>95%
buffer:	10 mM Tris-HCl, pH 8.0
preservative:	0.05% NaN ₃
isotype:	IgG ₁
concentration:	[lot specific] mg/mL
specificity:	Human CRP
cross-reactivity:	No cross reaction with ALB, TF, HGB, BIL and LDL in serum/plasma
pairing:	This antibody can be used in sandwich ELISA as the capture antibody in conjunction with A58016 or A58018 (as the detection antibody).
notes/applications:	ELISA or Lateral flow
storage:	2-8°C short term (<6 months); -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.