



## PRODUCT SPECIFICATION

**CATALOG NUMBER: A58060**

PRODUCT:	C-reactive protein (CRP) Monoclonal Antibody
IMMUNOGEN:	Human CRP
FORM:	DEAE- chromatography purified IgG
BUFFER:	0.015 M potassium phosphate buffer, 0.15M NaCl, pH 7.2
PRESERVATIVE:	0.1% NaN <sub>3</sub>
SUBCLASS:	IgG2a
STORAGE CONDITIONS:	2-8°C short term; -20°C long term storage
CROSS REACTIVITY:	None to other serum proteins
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.