



PRODUCT SPECIFICATION

CATALOG NUMBER: A58100

PRODUCT:	C-reactive protein (CRP) Monoclonal Antibody
IMMUNOGEN:	Human CRP
PURIFICATION:	Protein A sepharose
BUFFER:	PBS, pH 7.4
PRESERVATIVE:	0.1% NaN ₃
SUBCLASS:	IgG ₁
STORAGE CONDITIONS:	2-8°C
SPECIFICITY:	Human CRP
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precaution. There is no known biohazard associated with this product.