



PRODUCT SPECIFICATION

CATALOG NUMBER: A58110

PRODUCT:	C-reactive protein (CRP) Monoclonal Antibody
IMMUNOGEN:	Human CRP
PURIFICATION:	Protein A sepharose
BUFFER:	PBS, pH 7.4
PRESERVATIVE:	0.1% NaN ₃
SUBCLASS:	IgG _{2a}
STORAGE CONDITIONS:	2-8°C
SPECIFICITY:	Human CRP
PAIRING INFORMATION:	Forms a suitable sandwich when paired with A58080, A58120, A58130.
HAZARD/BIOHAZARD:	This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.