



PRODUCT SPECIFICATION

CATALOG NUMBER: A58120

| | |
|----------------------|--|
| PRODUCT: | C-reactive protein (CRP) Monoclonal Antibody |
| IMMUNOGEN: | Human CRP |
| PURIFICATION: | Protein A sepharose. |
| BUFFER: | PBS, pH 7.4 |
| PRESERVATIVE: | 0.1% NaN ₃ |
| SUBCLASS: | IgG _{2a} |
| STORAGE CONDITIONS: | 2-8°C |
| SPECIFICITY: | Human CRP |
| PAIRING INFORMATION: | Forms a suitable sandwich when paired with A58080, A58110, A58130. |
| HAZARD/BIOHAZARD: | This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product. |