



## PRODUCT SPECIFICATION

**CATALOG NUMBER: A58190**

PRODUCT:	Transferrin Monoclonal Antibody
IMMUNOGEN:	Human Transferrin
PURIFICATION:	Purified by chromatography on protein G Sepharose.
BUFFER:	PBS, Ph 7.4
PRESERVATIVE:	0.1% sodium azide
SUBCLASS:	IgG <sub>1</sub>
STORAGE CONDITIONS:	2-8°C
REACTIVITY:	Does not cross react with mouse, rabbit, or bovine transferrin. Does not cross react with human lactoferrin.
APPLICATION:	Detection of transferrin. Clone works in ELISA.
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.