



PRODUCT SPECIFICATION

CATALOG NUMBER: A58210

| | |
|---------------------|--|
| PRODUCT: | Human Transferrin Monoclonal Antibody |
| FORM: | Protein A purified in PBS buffer with pH 7.0 |
| PRESERVATIVE: | 0.02% NaN ₃ |
| SUBCLASS: | IgG ₁ |
| STORAGE CONDITIONS: | 2-8°C |
| HAZARD/BIOHAZARD: | This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product. |