



**PRODUCT SPECIFICATION**

**CATALOG NUMBER: A58220**

PRODUCT:	Human Transferrin Monoclonal Antibody
FORM:	DEAE-chromatography purified IgG
BUFFER:	0.015 M potassium phosphate buffer, 0.15 M NaCl, pH 7.2
PRESERVATIVE:	0.1% NaN <sub>3</sub>
SUBCLASS:	IgG <sub>1</sub>
STORAGE CONDITIONS:	2-8°C short term -20°C long term
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.