



PRODUCT SPECIFICATION

| | |
|----------------------|---|
| catalog number: | A58340 |
| product description: | Serum Amyloid A (SAA) Monoclonal Antibody, Purified |
| lot number: | [lot specific] |
| immunogen: | Purified Human Native SAA Protein |
| form: | Liquid |
| purification method: | Affinity chromatograph purified Fab; 0.2µm Filtered |
| purity: | >98% by SDS-PAGE |
| buffer: | PBS, pH 7.2 |
| preservative: | 0.05% NaN ₃ |
| concentration: | [lot specific] mg/mL |
| pairing: | For IT Assay use A58340/A53841/A58342 in 2/1/1 ratio and conjugate to latex For LF use A58340 as conjugate with A58341 and A58342 as capture |
| storage: | 2-8°C |
| hazard/biohazard: | This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product. |