



PRODUCT SPECIFICATION

catalog number:	A58342
product description:	Serum Amyloid A (SAA) Monoclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Purified Human Native SAA Protein
form:	Liquid
purification method:	Affinity chromatograph purified Fab; 0.2µm Filtered
purity:	>98% by SDS-PAGE
buffer:	PBS, pH 7.2
preservative:	0.05% NaN ₃
concentration:	[lot specific] mg/mL
pairing:	For IT Assay use A58340/A53841/A58342 in 2/1/1 ratio and conjugate to latex For LF: Use A58342 as label and pair with A58341 as capture
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.