



## PRODUCT SPECIFICATION

catalog number:	<b>A58412535P</b>
product description:	Phospho Tau 217 Monoclonal Antibody
lot number:	[lot specific]
source:	Murine
form:	Liquid
purity:	>95% (SDS PAGE)
buffer:	Phosphate Buffered Saline (PBS), pH 7.4
preservative:	0.1% NaN <sub>3</sub>
isotype:	IgG1
concentration:	[lot specific] mg/mL
expiration date:	[lot specific]
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.