



PRODUCT SPECIFICATION

CATALOG NUMBER: A60010

PRODUCT:	Influenza A Monoclonal Antibody
FORM:	Protein A purified Mouse MoAb
BUFFER:	10 mM sodium phosphate, 0.85% (w/v) NaCl, pH 7.2
PRESERVATIVE:	0.05% (w/v) sodium azide
FILTER:	0.2 um filtered solution
AFFINITY CONSTANT:	N/A
SUBCLASS:	IgG2a Kappa
STORAGE CONDITIONS:	2-8°C
SPECIFICITY:	Reacts with Influenza A subtypes H1N1, H2N2, H3N2. Reacts with nucleoprotein, shows no HI activity.
REFERENCE ASSAY:	Indirect fluorescent assay
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.