



PRODUCT SPECIFICATION

catalog number:	A60011
product description:	Influenza A Monoclonal Antibody, Purified
lot number:	[lot specific]
source:	Murine
form:	Liquid
purification method:	Protein A chromatography
purity:	>90%
buffer:	0.01M PBS, pH 7.2
preservative:	0.1% NaN ₃
isotype:	IgG _{2a}
concentration:	[lot specific] mg/mL
specificity:	This antibody is specific for nuclear protein of multiple strains, >50 different isolates from H1-H14.
pairing:	This antibody pairs well as the capture antibody with A60014 as the conjugate. This antibody also self pairs in lateral flow assays.
storage:	2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.