



PRODUCT SPECIFICATION

catalog number:	A60065
product description:	Influenza A Monoclonal Antibody
lot number:	[lot specific]
source:	Murine
form:	Liquid
purity:	≥95% by SDS-PAGE
buffer:	Citrate buffer, pH 6.0
preservative:	0.09% NaN ₃
isotype:	IgG ₁
concentration:	[lot specific] mg/mL
specificity:	Specific for Influenza A protein detection.
cross-reactivity:	No cross reactivity with Adenovirus or RSV.
pairing:	This antibody self-pairs.
expiration date:	[lot specific]
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.