



PRODUCT SPECIFICATION

catalog number:	A60068
product description:	Influenza A Monoclonal Antibody
lot number:	[lot specific]
source:	Murine ascites
form:	Liquid
purification method:	Protein A chromatography
purity:	>90%
buffer:	50 nM NaCl, 20 mM PBS, pH 7.4
preservative:	0.1% NaN ₃
isotype:	IgG _{2a}
concentration:	[lot specific] mg/mL
specificity:	Specific for Influenza A virus. Recognizes the nucleoprotein of numerous strains.
cross-reactivity:	No cross-reactivity with Parainfluenza, RSV, Adenovirus.
pairing:	This antibody can be used as labeling with A60069574P as coating antibody.
storage:	2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.