



PRODUCT SPECIFICATION

CATALOG NUMBER: A60120

PRODUCT:	Parainfluenza Type I Monoclonal Antibody
IMMUNOGEN:	Parainfluenza I C-35
FORM:	Mouse Ascites
SUBCLASS:	IgG _{2a} Kappa
STORAGE CONDITIONS:	4 C short term -20 C long-term. Repeated freeze-thaw will cause degradation of reagent.
SPECIFICITY:	Specific for the F0/F1 fusion protein common to all Parainfluenza I strains. No cross-reactivity to Parainfluenza II or III.